

MicroCap™ Depth Filter Capsules

Uniquely Flexible to Meet Your Processing Needs





Efficient, Cost-Effective Batch Processing

A Capsule Suite Customized to Fit Your Processing Needs

Eliminate Batch Pooling

Tired of expensive and inefficient pooling of multiple batches to sustain processing?

MicroCap Capsules offer the most flexible solution and the most consistent results for all of your batch processing needs.

MicroCap single-use capsules are a uniquely flexible line of disposable depth filter products designed to optimize scale-up and scale-down studies.

The MicroCap Capsule Suite with six different capsule sizes can meet your process volume requirements without requiring the pooling of multiple batches.

Lab-scale through clinical scale process volumes can be easily managed within the framework of cost-effective, efficient processing. Scale-up and scale-down studies can be efficiently managed through capsules sized and aligned to standard process volumes.

Advantages

- **Uniquely flexible** — options in capsule size, effective filter area and connection styles to suit your needs
- **Linear Scalability** — assurance of application and throughput from lab to production scale
- **Low hold-up volume** — reduced post-use rinsing volumes for product recovery
- **Completely disposable** — no cleaning or cleaning validation

Processes

- Mammalian cell cultures
- Bacteria, yeast, and insect cell lysates
- Vaccines
- Blood plasma proteins and serum
- Media

Applications

- Primary separations/prefiltration
- Secondary clarification
- Cell culture harvest
- Cell culture clarification
- DNA reduction
- Endotoxin reduction
- Host Cell Protein (HCP) reduction
- Protein aggregate removal
- Decolorization

Advanced Performance

An extensive range of enhanced depth filter medias has been developed by ErtelAlsop to meet the stringent requirements of the Biopharmaceutical industry and its unique filtration needs. Supported by comprehensive validation, ErtelAlsop's range of depth filter media enables reliable and efficient performance.

Reliability

All MicroCap depth filter capsules provide performance consistency and lot-to-lot traceability in easy-to-use formats. All MicroCap capsules are batch tested in order to meet all quality requirements and meet all applicable USP requirements including Class VI Biological Testing for Plastics.

Scalability

The versatility of MicroCap capsules enhances filtration efficiencies of laboratory, pilot and small-scale processes. The range of capsules and respective filter area enable both linear scalability in performance and as well as through the range of capsule sizes.




Filter Area

MicroCap Capsule	Single Layer (cm ²)*	Double Layer (cm ²)**
MC1	26	26
MC2	170	80
MC5	330	170
MC10	960	500
MC20	1920	1000
MC30	2880	1500

* Single-layer media such as MicroMedia®, XL Series™, and MicroClear™
 ** Double-layer media such as DXL Series

Approximate Batch Size



Efficient, Cost-Effective Processing

Available in a variety of sizes and configurations, MicroCap Capsules enable single-batch processing without costly pooling of multiple batches. This variety in capsule sizes and configurations also provides adequate filter area to complete entire processing runs in an efficient and cost effective manner.

Consistent Processing

MicroCap Capsules offer processing consistency through manufacturing controls of the depth filter media, creating the lot-to-lot consistency of the filter media required for performance reliability. These manufacturing controls, coupled with capsule design that utilizes the same internal structure and flow path, enables consistent processing from run to run and from capsule to capsule.

Predictable Processing

Consistent product design, internal structure and flow path coupled with reliable filter media performance enable predictive scalability throughout the range in capsule sizes.





Capsules for Process Development and Lab Scale Processing

MicroCap capsules are available in 6 different sizes (1 in., 2.5 in., 5 in., 10 in., 20 in. and 30 in. in height) ranging in effective filter area from 23 cm² to 2880 cm² with single-layer media (23 cm² to 1500 cm² with double-layer media). This enables processing without the necessity of pooling multiple batches of material to efficiently utilize the appropriate single-use capsule.

The MC1 capsule* is designed to suit your testing needs as a screening tool for depth filter media selection. The MC1 provides optimal throughput and predictive, scalable results from small batch sizes ranging from approximately 0.5 liters to 3 liters.

Positioned for use during laboratory development and optimization of scale-up and scale-down studies, MC1 capsules, with 23 cm² of effective filter area (EFA), are effective tools for quickly determining the correct depth filter media to suit your processing needs. With the same filter media, flow paths and design as the larger capsules in this series, the MC1 capsule allows for simplified process development studies.

The MC2 and MC5 capsules are intended for use in processes batch sizes typical of early stage laboratory to process development ranging from approximately 3 to 15 liters. Available with a large variety of depth filter media options and optional inlet/outlet configurations, MicroCap capsules provide simple, reliable capsule performance suited to your process needs.

*The MC1 is also available in the MicroCap Laboratory Cabinet (pg. 11), or reference: Technical Bulletin Cab-13.

MC1



MC2



MC5





Capsules for Lab and Production Scale Processing

Intended for larger volume processing (10 to 50 or more liters), the MC10, MC20, and MC30 MicroCap capsules share the same internal structure, flow path, and robust design as the smaller capsules in the MicroCap line. These T-style capsules are easily manifolded to seamlessly integrate additional filter area or establish a fully encapsulated filter train. With effective filter areas of 960 cm², 1920 cm² and 2880 cm² in single-layer format (500 cm², 1000 cm² and 1500 cm² in double-layer format) the MC10, MC20 and MC30 MicroCap capsules are ideally sized for larger process development though production scale needs.

MC10



MC20




MC30



Common Flow Path for Predictable Scaleup

MicroCap Capsules use the same internal cartridge design for processing consistency. Their consistent structure and flow path provide enhanced and predictable linear scale-up for all capsule sizes.





MicroCap Filter Media

MicroMedia

Pharmaceutical grade MicroMedia Filter Sheets contain a high performance filter media consisting of a balance of cellulose fibers, diatomaceous earth (DE) filter aid, and a wet strength resin. MicroMedia Filter Sheets provide consistent performance across a variety of applications.

XL Series

ErtelAlsop's MicroMedia XL Series depth filter media incorporates a high performance grade of diatomaceous earth (DE) for greatly enhanced filter performance. In addition, its high purity, low extractables and rigorous quality control make it the leading product for depth filtration in critical applications such as pharmaceuticals and biologics. The reduced levels of metals/contaminants result in increased permeability compared to conventional DE contribute to its purity and usability as a filter for products that must meet the highest standards.

DXL Series

Consisting of two distinct layers, the double-layer configuration of DXL Series depth filter media maximizes contaminant loading within the body of the filter structure for prolonged filter life. DXL Series depth filter media efficiently manage process streams containing high solids that are often associated with low viability cell cultures. Combining any two layers of ErtelAlsop's XL Series media allows the design of a filter solution to address the unique characteristics of a process stream and to improve process optimization.

MicroClear

Consisting of immobilized activated carbon and cellulose fibers, ErtelAlsop's MicroClear depth filter media provides superior adsorption characteristics over traditional stirred tank or packed bed carbon applications. Designed for use in a single-pass process, the large surface area of carbon to the process stream enables highly efficient color or odor removal. Unlike other depth filter media, designed for mechanical particle removal, MicroClear depth filter media are designed for adsorption of colors, odors or other soluble contaminants.

MicroClear Carbon per Unit Area

Effective Filter Area	Amount of Carbon (grams)
26 cm ²	2
170 cm ²	15
330 cm ²	29
960 cm ²	83
1920 cm ²	166
2880 cm ²	249


MicroClear Activation Methods and Typical Applications

Media Grade	Activation Method	Typical Application
MC55P	Steam Activated	General Applications
MC55CP	Chemically Activated	Decolorization in Pharmaceutical Applications
MC55GP	Steam Activated	Fine Chemical and Pharmaceutical intermediates






Fine to Medium Filter Media Grades




MicroMedia or XL Series filter sheets for the clarification of very fine to medium (0.25 - 2.0 micron) sized contaminant from process streams. Available in all six sizes of MicroCap single-use capsules.

Medium to Coarse Filter Media Grades



Process streams containing medium to coarse (1.0 - 25 micron) sized contaminant, XL Series depth filter media grades provide enhanced throughput and lower extratables. Available in all MicroCap single-use capsules.


Double Layer (DXL) Filter Media Grades




Developed for process streams with high solids removal needs, large particle distributions, or reduction of staged filtration. The DXL grades combine two distinct layers of XL series depth filter media in a single capsule. Scalable results in all six MicroCap single-use capsule sizes.




MC1




MC2




MC5




MicroCap Dimensions (mm)						
Inlet/Outlet Fittings	MC1		MC2		MC5	
	A	B	A	B	A	B
Hose Barb	-	79	80	148	137	205
Sanitary Clamp	-	-	80	126	137	183



MC10



MC20



MicroCap Dimensions (mm)									
Inlet/Outlet Fittings	MC10			MC20			MC30		
	A	B	C	A	B	C	A	B	C
Sanitary Clamp	304	356	127	547	599	127	790	842	127





Technical Specifications

Materials of Construction	
Capsule Shell and internals	Polypropylene
Capsule Filter Media	Cellulose fibers Diatomaceous earth Powdered activated carbon* Resin
O-rings	Silicone
Maximum Operating Pressure	
MicroCap MC1	2.5 barg (36 psig) @ 22 °C
All Other MicroCap Capsules	5.5 barg (80 psig) @ 22 °C
Maximum Differential Pressure	
MicroCap MC1	2 bard (30 psid)
All Other MicroCap Capsules	2.4 bard (35 psid)
Sterilization - Autoclave	
MicroCap MC1	1 cycle @ 121 °C for 30 minutes
All Other MicroCap Capsules	2 cycles @ 125 °C for 60 minutes

* Carbon with MicroClear Capsules Only

Hold-up Volume

MicroCap Capsules	Total Void Volume (mL)*	Holdup Volume (mL)**
MC1	-	-
MC2	176	68
MC5	299	142
MC10	1080	392
MC20	2160	784
MC30	3240	1176

* **Total Void Volume** = Total upstream and downstream void volume with internal cartridge installed

** **Holdup Volume** = Liquid remaining in depth filter media post blow down.



MicroCap Laboratory Cabinets

With single-use technology on the rise, our new MicroCap Laboratory Cabinet has been created to simplify filtration and separation trials for research and development, as well as, process development laboratories. This attractive and well organized Cabinet easily stores the filters you require.

There are five cabinet options to choose from each containing MicroCap MC1 capsules with 23 cm² (3.48 in.2) of effective filtration area.

Whether your needs are clarification, cell harvest or color removal there is a cabinet designed to meet your needs. The easy to use capsules allow users to test a comprehensive range of depth filter media and efficiently determine the right solution for their process needs. For more information, please reference: Technical Bulletin Cab-13.



MicroCap Media Grades

Media Series	Media Grades	Nominal Rating (µm)	Format	Filter Properties/Material	
MicroMedia	M954P	0.25 – 1.0	Single Layer	Cellulose with filter aid	
	M854P	0.3 – 1.25			
	M704P	0.45 – 1.5			
XL	M953P	0.25– 1.0		Cellulose with high purity filter aid	
	M853P	0.3 – 1.25			
	M703P	0.45 – 1.5			
	M503TP	0.8 – 2.75			
	M503P	1.0 – 3.0			
	M453P	2.5– 6.0			
	M403P	5.0 – 12.0			
DXL*	M103P	10.0 – 17.0		Double Layer	Cellulose with high purity filter aid
	M053P	15.0 – 20.0			
	B9E9	0.25 – 1.0			
	B6E9	0.25 – 2.75			
	B5E8	0.3 – 3.0			
	B4E7	0.45 – 6.0			
MicroClear	B2E6	0.8 – 17.0	Single Layer	Cellulose with activated carbon	
	B1E5	1.0 – 20.0			
	MC55P	Steam Activated			
	MC55CP	Chemical Activated			
	MC55GP	Steam Activated			

* Other media combinations available upon request



MicroCapTM

Single-Use Capsules

MicroCap Ordering Information

Example: M053PCAP05MTB is a MC5 capsule with M053P media, 1/2 in. sanitary clamp inlet/outlet, and vent valves.

Media Grade	Type	Size	Inlet/Outlet	Description	Vent	Description
See Media Grades Table	CAP	Z1	H	Stepped Hose Barb	O	Luer Lock
		02	3H	3/8 in. Hose Barb	B	Vent Valve
		05	MT	1/2 in. Sanitary Clamp		
		10	TC	1 1/2 in. Sanitary Clamp	C	Inlet/Outlet 1/4 in. Bleed Valve
		20				
		30				

MicroCap Laboratory Cabinet Ordering Information

Part #	Media Types Included	Number of Capsules Held	Material of Construction	Cabinet Dimensions
MCSL01	MicroMedia, MicroMedia: XL Series	48	Polypropylene	10.2"w x 14.5"d x 12.8"h
MCSL02	MicroClear	15	Polypropylene	10.2"w x 14.5"d x 4.5"h
MCSL03	MicroMedia, MicroMedia: XL Series, MicroClear	63	Polypropylene	10.2"w x 14.5"d x 16.8"h
MCDL01	MicroMedia: DXL Series	48	Polypropylene	10.2"w x 14.5"d x 12.8"h
MCCC01	Customizable up to 6 media grades	48	Polypropylene	10.2"w x 14.5"d x 12.8"h



ENFISO GmbH
Weinbergstrasse 4
8447 Dachsen
+41 79 752 50 82
info@enfiso.ch | www.enfiso.ch

Because of developments related to products, systems and services, the products, data and procedures are subject to change without notice.

Please consult your ErtelAlsop representative or visit www.ErtelAlsop.com to verify that this information remains valid.

MicroCap, MicroMedia, and MicroClear are trademarks of ErtelAlsop.

© Copyright 2013, ErtelAlsop

Technical Bulletin MCI-13 Rev C

